

## Sample Packing Tips

---

### Tube Shipping Preparation

Tubes should be secured with plastic wrap prior to shipping. This keeps tubes in order and greatly reduces the chance of tubes breaking during transit.



Example of tubes secured in plastic wrap.

Example of unsecured tubes upon their arrival in the lab.  
Tubes may arrive broken or damaged otherwise.



Armor Animal Health's Lab is located 6040 N Cutter Cir # 317, Portland, OR 97217

For more information call Toll-Free: 877.278.1344

Copyright © 2021 Armor Animal Health. All rights reserved.