

# Stocker/Backgrounder BVD PI Screening Protocol



The screening protocol is designed to keep Bovine Viral Diarrhea (BVD) persistently infected (PI) calves from causing production losses for the stocker-backgrounder. Identifying and eliminating PIs early can realize benefits of reduced disease and lowered treatment costs for the operation. Selling backgrounded cattle as tested BVD PI negative increases their value to the feed yard buyer. All beef operations can benefit from a BVD control plan involving testing, proper vaccination, biosecurity and good animal handling and welfare management.

1. All cattle should have an ear notch or blood sample taken on arrival.
2. All cattle need to be identified individually and matched to the sample taken.
3. Samples should be sent frequently to the lab for analysis. If possible, keep the sampled group in a separate pen without contact with other animals or shared water tanks and feed sources until the test results are known.
4. Positive PI calves should be removed immediately from the group as PIs shed virus in large amounts continuously and spread the virus to other calves.
5. PI positive calves should be isolated in a separate part of the yard. Care should be taken to avoid spreading the virus from these calves to other calves in the yard.
6. A confirmatory sample should be taken from the PI positive calf 3 weeks after the original test.
7. Positively confirmed PI cattle can be railed after processing and treatment withdrawals are clear. Some operations elect to finish PI animals separately. Approximately 30-50% of PIs will die prior to reaching optimal finish weight. We do not recommend marketing the animal.
8. Negative second test calves were acutely infected and can be returned to the home pen or group. An acute calf is one that was infected by a PI, but was able to clear the infection.
9. A system should be developed to identify which buyers or sale barns are sending positive calves. Information should be sent to them about ways that they can implement a BVD control program. Armor Animal Health is committed to working with the industry to assist in a program designed to stop PIs from limiting production.

Armor Animal Health's Lab is located 6040 N Cutter Cir # 317, Portland, OR 97217

For more information call Toll-Free: 877.278.1344